
CoTACIT step by step Protocol

From the He lab @PKU

Ver.2025, Min Liu

Buffers:

Reagents	Composition
Wash Buffer	1 mL 1 M HEPES-KOH pH 7.5, 1.5 mL 5 M NaCl, 12.5 μ L 2 M spermidine, 500 μ L sodium butyrate, 100 μ L 5% digitonin, 125 μ L 20% Triton X-100, 1 \times cocktail, and the final volume to 50 mL with ddH ₂ O
High-salt Wash Buffer	1 mL 1 M HEPES-KOH pH 7.5, 3 mL 5 M NaCl, 12.5 μ L 2 M spermidine, 500 μ L sodium butyrate, 100 μ L 5% digitonin, 125 μ L 20% Triton X-100, 1 \times cocktail, and the final volume to 50 mL with ddH ₂ O
Antibody Buffer	Mix 4 μ L 0.5 M EDTA with 1 mL Wash-buffer
Tagmentation Buffer	20 mM TAPS-NaOH pH 8.3, 20 mM MgCl ₂ , 1 \times cocktail, 20 mM sodium butyrate, and 2 mM PMSF
Lysis Buffer	10 mM Tris-HCl pH 8.5, 0.05% SDS and 0.1 mg/mL proteinase K

Notes: The hyperactivity of protein A-Tn5 (PAT) is most critical in this protocol to achieve extremely high genome coverage in single cells. Most commercially available PAT products are not so good in our hands. He lab is happy to provide a reasonable amount of in-house PAT upon request.

1. Take out methanol-fixed cells from -80°C freezer and put on ice for at least 15 minutes for rehydration.
2. Wash cells twice with 180 μ L Wash Buffer (or 0.1% BSA/PBS) to remove residual methanol. For each time, cells were collected at an appropriate speed for 3 minutes at 4°C and 10-20 μ L Wash Buffer was retained to reduce cell loss.

First round barcoding:

3. Resuspend cells with 100 μ L Antibody Buffer containing 0.5 μ g primary antibodies

(for example, anti-H3K27ac antibody).

4. Put the tube on a rotator and incubate at 4°C for 4 hours.
5. Collect cells at an appropriate speed for 3 minutes at 4°C.
6. Resuspend cells with 180 µL Wash Buffer and rotate cells at 4°C for 5 minutes to wash out free antibodies.
7. Collect cells at an appropriate speed for 3 minutes at 4°C.
8. Resuspend cells with 100 µL High-salt Wash Buffer containing 3 µg/ml PAT-T5-1/T7-1 enzyme.
9. Put the tube on a rotator and incubate at 4°C for 60 minutes.
10. Wash cells twice with 180 µL High-salt Wash Buffer and rotate cells at 4°C for 5 minutes to wash out free PAT enzyme.
11. Collect cells at an appropriate speed for 3 minutes at 4°C. Aspirate the supernatant and retain about 10 µL High-salt Wash Buffer.
12. Resuspend cells with 10 µL Tagmentation Buffer and incubate at 37°C for 60 minutes.
13. Add 20 µL 40 mM EDTA and incubate at room temperature for 15 minutes to stop the reaction.
14. Wash cells three times with 0.1% BSA/PBS (or Wash Buffer) and rotate cells 4°C for 5 minutes to wash out free PAT enzyme.

Second round barcoding:

15. Resuspend cells with 100 µL Antibody Buffer containing 0.5 µg primary antibodies (for example, anti-H3K27me3 antibody).
16. Put the tube on a rotator and incubate at 4°C for 4 hours.
17. Collect cells at an appropriate speed for 3 minutes at 4°C.
18. Resuspend cells with 180 µL Wash Buffer and rotate cells at 4°C for 5 minutes to wash out free antibodies.
19. Collect cells at an appropriate speed for 3 minutes at 4°C.
20. Resuspend cells with 100 µL High-salt Wash Buffer containing 3 µg/ml PAT-T5-2/T7-2 enzyme.
21. Put the tube on a rotator and incubate at 4°C for 60 minutes.
22. Wash cells twice with 180 µL High-salt Wash Buffer and rotate cells at 4°C for 5 minutes to wash out free PAT enzyme.
23. Collect cells at an appropriate speed for 3 minutes at 4°C. Aspirate the supernatant and retain about 10 µL High-salt Wash Buffer.
24. Resuspend cells with 10 µL Tagmentation Buffer and incubate at 37°C for 60

minutes.

25. Add 20 μ L 40 mM EDTA and incubate at room temperature for 15 minutes to stop the reaction.
26. Wash cells three times with 0.1% BSA/PBS (or Wash Buffer) and rotate cells 4°C for 5 minutes to wash out free PAT enzyme.

Third round barcoding:

27. Resuspend cells with 100 μ L Antibody Buffer containing 0.5 μ g primary antibodies (for example, anti-H3K9me3 antibody).
28. Put the tube on a rotator and incubate at 4°C for 4 hours.
29. Collect cells at an appropriate speed for 3 minutes at 4°C.
30. Resuspend cells with 180 μ L Wash Buffer and rotate cells at 4°C for 5 minutes to wash out free antibodies.
31. Collect cells at an appropriate speed for 3 minutes at 4°C.
32. Resuspend cells with 100 μ L High-salt Wash Buffer containing 3 μ g/ml PAT-T5-3/T7-3 enzyme.
33. Put the tube on a rotator and incubate at 4°C for 60 minutes.
34. Wash cells twice with 180 μ L High-salt Wash Buffer and rotate cells at 4°C for 5 minutes to wash out free PAT enzyme.
35. Collect cells at an appropriate speed for 3 minutes at 4°C. Aspirate the supernatant and retain about 10 μ L High-salt Wash Buffer.
36. Resuspend cells with 10 μ L Tagmentation Buffer and incubate at 37°C for 60 minutes.
37. Add 20 μ L 40 mM EDTA and incubate at room temperature for 15 minutes to stop the reaction.
38. Wash cells three times with 0.1% BSA/PBS (or Wash Buffer) and rotate cells 4°C for 5 minutes to wash out free PAT enzyme.

Library construction:

39. Pre-rinse the 96-well plates with 0.1% BSA/PBS to avoid loss of DNA fragments.
40. And add 2 μ L Lysis Buffer to each well of the pre-rinsed plates.
41. Pick single cell into a well of a 96-well plate with a mouth pipette under a microscope.
42. Add 5 μ L mineral oil (Sigma) to each well and incubate the plates at 55°C for 15 minutes to release DNA fragments.
43. Add 0.5 μ L of 10 mM PMSF and 1 μ L of 0.9% Triton X-100 to each well to deactivate protease K and quench SDS.

-
44. Add 17 μL PCR mix (0.2 μL KAPA HiFi HotStart DNA polymerase, 4 μL 5 \times KAPA High-GC buffer, 0.5 μL 10 mM dNTP mix, 0.5 μL 25 mM MgCl_2 , 11.3 μL H_2O and 0.5 μL 50 μM first-round primer mix) to each well.
 45. Perform PCR with the following program: 1 cycle of 72°C for 5 minutes; 1 cycle of 95°C for 3 minutes; 17 cycles of 98°C for 20 seconds, 65°C for 30 seconds, 72°C for 1 minutes; 1 cycle of 72°C for 5 minutes; and hold at 4°C.
 46. Add 0.25 μL ExoI (NEB) enzyme and incubate plates at 37°C for 60 minutes to digest excess primers.
 47. Incubate plates at 72°C for 20 minutes to deactivate ExoI enzyme.
 48. 10 μL second-round PCR mix (0.1 μL KAPA HiFi HotStart DNA polymerase, 2 μL 5 \times KAPA High-GC buffer, 0.25 μL 25 mM MgCl_2 and 6.4 μL ddH₂O) was added to each well containing 0.5 μL 10 mM Truseq i5 and 0.5 μL Truseq i7 index primers.
 49. Perform PCR with the following program: 1 cycle of 95°C for 3 minutes; 5 cycles of 98°C for 20 seconds, 65°C for 30 seconds, 72°C for 1 minutes; 1 cycle of 72°C for 5 minutes; and hold at 4°C.
 50. Purify the DNA fragments with 1 \times AMPure XP beads (Beckman) once, and 200-1,000 bp fragments were selected with 0.5 \times + 0.5 \times AMPure XP beads.
 51. The libraries were sequenced with paired-end 150-bp reads on a NovaSeq 6000 platform (Illumina) or BGI T7 or MGI2000 platform.

Table 1. The sequences of PAT-T5/T7 used in CoTACIT.

Primer	Sequence
T5-1	TCGTCGGCAGCGTCTCCACGCTATAGCCTGCGATCGAGGACGGCAGAT GTGTATAAGAGACAG
T5-2	TCGTCGGCAGCGTCTCCACGCATAGAGGCGCGATCGAGGACGGCAGA TGTGTATAAGAGACAG
T5-3	TCGTCGGCAGCGTCTCCACGCCCTATCCTGCGATCGAGGACGGCAGAT GTGTATAAGAGACAG
T7-1	GTCTCGTGGGCTCGGCTGTCCCTGTCCCGAGTAATCACCGTCTCCGCCT CAGATGTGTATAAGAGACAG
T7-2	GTCTCGTGGGCTCGGCTGTCCCTGTCCCTCTCCGGACACCGTCTCCGCCT CAGATGTGTATAAGAGACAG
T7-3	GTCTCGTGGGCTCGGCTGTCCCTGTCCAATGAGCGCACCGTCTCCGCCT CAGATGTGTATAAGAGACAG
common annealing primer	5Phos/CTGTCTCTTATACACATCT

Table 2. The sequences of Truseq i5 and Truseq i7 index primers used in CoTACIT.

Primer	Sequence
Truseq-501	AATGATACGGCGACCACCGAGATCTACACT TATAGCCT ACACTCTTTCCCTACACGA CGCTCTTCCGATCT
Truseq-502	AATGATACGGCGACCACCGAGATCTACAC ATAGAGGC ACACTCTTTCCCTACACGA CGCTCTTCCGATCT
Truseq-503	AATGATACGGCGACCACCGAGATCTACAC CCTATCCT ACACTCTTTCCCTACACGAC GCTCTTCCGATCT
Truseq-504	AATGATACGGCGACCACCGAGATCTACAC GGCTCTGA ACACTCTTTCCCTACACGA CGCTCTTCCGATCT
Truseq-505	AATGATACGGCGACCACCGAGATCTACAC AGGCGAAG ACACTCTTTCCCTACACGA CGCTCTTCCGATCT
Truseq-506	AATGATACGGCGACCACCGAGATCTACAC TAATCTTA ACACTCTTTCCCTACACGA CGCTCTTCCGATCT
Truseq-507	AATGATACGGCGACCACCGAGATCTACAC CAGGACGT ACACTCTTTCCCTACACGA CGCTCTTCCGATCT
Truseq-508	AATGATACGGCGACCACCGAGATCTACAC GTACTGAC ACACTCTTTCCCTACACGA CGCTCTTCCGATCT
Truseq-509	AATGATACGGCGACCACCGAGATCTACACT TGCTTG ACACTCTTTCCCTACACGAC GCTCTTCCGATCT
Truseq-510	AATGATACGGCGACCACCGAGATCTACAC GAGAGGTT ACACTCTTTCCCTACACGA CGCTCTTCCGATCT
Truseq-511	AATGATACGGCGACCACCGAGATCTACAC ACCTGGTT ACACTCTTTCCCTACACGA CGCTCTTCCGATCT
Truseq-512	AATGATACGGCGACCACCGAGATCTACAC AAGCGGAA ACACTCTTTCCCTACACGA CGCTCTTCCGATCT
Truseq-513	AATGATACGGCGACCACCGAGATCTACAC CGGAACAA ACACTCTTTCCCTACACGA CGCTCTTCCGATCT
Truseq-514	AATGATACGGCGACCACCGAGATCTACAC GGTAAGCT ACACTCTTTCCCTACACGA CGCTCTTCCGATCT
Truseq-515	AATGATACGGCGACCACCGAGATCTACACT TGTGGCAT ACACTCTTTCCCTACACGA CGCTCTTCCGATCT

Primer	Sequence
Truseq-516	AATGATACGGCGACCACCGAGATCTACACACTACGGAACACTCTTTCCTACACGA CGCTCTTCCGATCT
Truseq-701	CAAGCAGAAGACGGCATACGAGATCGAGTAATGTGACTGGAGTTCAGACGTGTGC TCTTCCGATCT
Truseq-702	CAAGCAGAAGACGGCATACGAGATTCTCCGGAAGTGACTGGAGTTCAGACGTGTGCT CTTCCGATCT
Truseq-703	CAAGCAGAAGACGGCATACGAGATAATGAGCGGTGACTGGAGTTCAGACGTGTGC TCTTCCGATCT
Truseq-704	CAAGCAGAAGACGGCATACGAGATGGAATCTCGTGACTGGAGTTCAGACGTGTGCT CTTCCGATCT
Truseq-705	CAAGCAGAAGACGGCATACGAGATTCTGAATGTGACTGGAGTTCAGACGTGTGCT CTTCCGATCT
Truseq-706	CAAGCAGAAGACGGCATACGAGATACGAATTCGTGACTGGAGTTCAGACGTGTGCT CTTCCGATCT
Truseq-707	CAAGCAGAAGACGGCATACGAGATAGCTTCAGGTGACTGGAGTTCAGACGTGTGCT CTTCCGATCT
Truseq-708	CAAGCAGAAGACGGCATACGAGATGCGCATTAAGTGACTGGAGTTCAGACGTGTGCT CTTCCGATCT
Truseq-709	CAAGCAGAAGACGGCATACGAGATCATAGCCGGTGACTGGAGTTCAGACGTGTGC TCTTCCGATCT
Truseq-710	CAAGCAGAAGACGGCATACGAGATTTCGCGGAAGTGACTGGAGTTCAGACGTGTGCT CTTCCGATCT
Truseq-711	CAAGCAGAAGACGGCATACGAGATGCGCGAGAAGTGACTGGAGTTCAGACGTGTGC TCTTCCGATCT
Truseq-712	CAAGCAGAAGACGGCATACGAGATCTATCGCTGTGACTGGAGTTCAGACGTGTGCT CTTCCGATCT
Truseq-713	CAAGCAGAAGACGGCATACGAGATAGGAGGAAAGTGACTGGAGTTCAGACGTGTGC TCTTCCGATCT
Truseq-714	CAAGCAGAAGACGGCATACGAGATAGCAAGCAAGTGACTGGAGTTCAGACGTGTGC TCTTCCGATCT
Truseq-715	CAAGCAGAAGACGGCATACGAGATTCATCACCGTGACTGGAGTTCAGACGTGTGCT CTTCCGATCT

Primer	Sequence
Truseq-716	CAAGCAGAAGACGGCATAACGAGATCGTAGGTTGTGACTGGAGTTCAGACGTGTGCT CTTCCGATCT
Truseq-717	CAAGCAGAAGACGGCATAACGAGATTCAGATCCGTGACTGGAGTTCAGACGTGTGCT CTTCCGATCT
Truseq-718	CAAGCAGAAGACGGCATAACGAGATCGTGATCAGTGACTGGAGTTCAGACGTGTGCT CTTCCGATCT
Truseq-719	CAAGCAGAAGACGGCATAACGAGATAGTCGTTGTGACTGGAGTTCAGACGTGTGCT CTTCCGATCT
Truseq-720	CAAGCAGAAGACGGCATAACGAGATGAACGCTTGTGACTGGAGTTCAGACGTGTGCT CTTCCGATCT
Truseq-721	CAAGCAGAAGACGGCATAACGAGATTACGCCTTGTGACTGGAGTTCAGACGTGTGCT CTTCCGATCT
Truseq-722	CAAGCAGAAGACGGCATAACGAGATCTATCAGGTGACTGGAGTTCAGACGTGTGCT CTTCCGATCT
Truseq-723	CAAGCAGAAGACGGCATAACGAGATTCTTCTGCGTGACTGGAGTTCAGACGTGTGCT CTTCCGATCT
Truseq-724	CAAGCAGAAGACGGCATAACGAGATGCTGGATTGTGACTGGAGTTCAGACGTGTGCT CTTCCGATCT